

Supplemental Table

Tab S1. Primers used for real-time PCR

| Genes | Primers |
|--------------------|--|
| NF- κ B p65 | Forward: 5'-ACTGCCAACAGCAGATGGCCCA Reverse: 5'-GCACCAGGTAGTCCACCATGGGAT |
| STAT3 | Forward: 5'-GGCCCCTCGTCATCAAGA Reverse: 5'-TTTGACCAGCAACCTGACTTTAGT |
| Bcl-2 | Forward: 5'-ATCGCCCTGTGGATGACTGAG Reverse: 5'-CAGCCAGGAGAAATCAAACAGAGG |
| GAPDH | Forward: 5'-AAGGACTCATGACCACAGTCCAT Reverse: 5'-CCATCACGCCACAGTTTCC |

Supplemental Figure Legend

Fig S1. Expression levels of NF- κ B p65, STAT3 and Bcl-2 in SKOV3 and SKOV3/DDP cells. SKOV3 and SKOV3/DDP cells were harvested followed by lysate preparation and Western blot analysis with antibodies indicated. GAPDH was used as a loading control.

