

Supplementary material

Figure 1. BIX-01294 induces apoptosis in human Molt-4 and Jurkat cells

- A. Molt-4 and Jurkat cells were incubated with the indicated concentrations of BIX-01294 for 36 h, cells were harvested and subsequently stained with annexin-V-PE and 7-AAD and analyzed by flow cytometry for apoptosis.
- B. Human PBMC were incubated with the indicated concentrations of BIX-01294 for 36 h, cells were harvested and subsequently stained with annexin-V-PE and 7-AAD and analyzed by flow cytometry for apoptosis.
- C. Molt-4 and Jurkat cells were incubated with the indicated concentrations of BIX-01294 for 48 h. Then whole cells were harvested and subjected to western blot using caspase 8, caspase 3, and PARP antibodies.

